
# Setting Up A Run on HiSeq Control Software (HCS2.0)

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\*Assumes High Output post-run maintenance wash has already been performed (~2 hours 40 minutes)

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A Welcome Screen Menu Button  
 B Interface Panel for Flow Cell A  
 C Interface Panel for Flow Cell B

3

## HCS Rapid Run Run Configuration: Wash

- ▶ Wash
  - Water またはMaintenance が選択可能
  - Water : MilliQのみ 15min程度
  - Maintenance: MilliQ → 1N NaOH → MilliQ 40min程度
  - Maintenance washの前にガasketの交換
  - 1日以上Idle状態のときwater washの実施が必須。

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

9

## Hiseq Rapid Run Run Configuration : Recipes

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10

# Hiseq Rapid Run Run Configuration : Sample sheet

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illumina®

11

# Hiseq Rapid Run Run Configuration : Reagents

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12

# Hiseq Rapid Run Run Configuration : Review

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13

# Hiseq Rapid Run Pre-Run Setup: Load Reagent

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14

# Hiseq Rapid Run Pre-Run Setup: Load Priming FC

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15

# Hiseq Rapid Run Initiate Run: Load Seq FC

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illumina®

17

# Hiseq Rapid Run Initiate Run: Close Door

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## FBR Guidelines\* (Rapid Run)

First Base Report (FBR) \ <run_folder&gt;\first_base_report< th=""><th>Expected value</th></run_folder&gt;\first_base_report<>	Expected value
Cluster density	400-600 K/mm2
C Intensity	>20000
A Intensity	>15000
G Intensity	>10000
T Intensity	>15000
Focus Score (Avg)	>55

\*Important: Values are only a guideline and each run should be assessed individually!

## FBR example

### First Base Report

120827\_D0003\_0438\_AH00GHADXX

### Top Surface

Metric	Lane 1	Lane 2
Cluster Density (k/mm2)	560.16	554.15
A Intensity	48517	47676
C Intensity	73616	72221.5
G Intensity	25786	26065.5
T Intensity	43358.5	43446.5
A Focus Score	82.3	82.19
C Focus Score	80.8	80.58
G Focus Score	82.94	83.11
T Focus Score	81.94	81.96

illumina

19

# Hiseq Rapid Run Initiate Run: Post Wash

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
illumina®

22

## Hiseq High Output Run Run Configuration : Flow Cell Set up

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illumina®

23

## Hiseq High Output Run Run Configuration : Advanced

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
24




# Hiseq High Output Run Run Configuration : Recipe

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25

# Hiseq High Output Run Run Configuration : Sample Sheet

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Please review the displayed values and select "NEXT" to continue.

**Flow Cell Setup**

Flow Cell ID: 8000000  
Equipment: ihi100000000  
User Name: ihiuser  
Flow Cell Type: HiSeq Flow Cell v3

**Storage**

Execution: ...  
MyName Username: N/A  
Output Folder: C:\Users\...  
Zip BCLs:   
Bin Q Scores:   
Save Aux Files: Save All Thumbnails

**Advanced**

Control Lane: None  
Confirm 1st Base:   
Keep Integrity Files: ...  
Align to Prox. Lanes: All Lanes

**Sample Sheet**

Sample Sheet: N/A

Navigation: BACK, CANCEL, HELP, NEXT

illumina®

28

# Hiseq High Output Run Pre-Run Setup: Load Reagent

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illumina®

29

# Hiseq High Output Run Pre-Run Setup: Load FC, Prime

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33

## Hiseq High Output Run Post wash: Maintenance wash

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